


# KT 625

## 6" (152mm) Temperature Chart Recorder

Our 6" (152mm) line of temperature chart recorders feature a sleek look and feel, precision pen movement, k-thermocouple probes and 8 user selectable temperature ranges to choose from. The KT625 has a digital display, two probes, relay contacts and alarms. The 6" (152mm) is a great size for monitoring refrigerators and freezers.



### Specification

- Alarm Type Audio/Visual on High/Low Temp Limits
- Ambient Operating RH Conditions 0-95% RH non-condensing
- Ambient Operating Temperature Conditions 32 to 122F (0 to 50C)
- Approvals RoHS, CE
- Battery Backup 9V Battery Back-up. Powers unit for up to 2 days.
- Battery Life (Avg) 48 Hours (Back-up Only)
- Calibration Options N100, N300, N400
- Channels CH1: Temperature
- Chart Rotation Speed 24 Hour or 7 Day
- Dimensions 3.25 x 7.56 x 7.56in
- Display Dimensions Height 1.18
- Display Dimensions Width 0.59
- Display Resolution 0.1 Degree from -9.9 to 199.9; 1 Degree resolution above 199 and below 9.9
- Display Type LCD
- Enclosure Rugged Black ABS Case, Polycarbonate Door, Aluminum Dial
- IP Rating 21
- Included Accessories AC Adapter, Red Pen, 9V Battery, Quick Start Guide and 4' Bead-Wire Probe. \*Charts sold separately.
- Keypad Functions Power On/Off, Pen Home, Alarm
- Mounting Options Free Standing or Wall Mounted
- Note on Relay Type Relay does not function while unit is running on only battery power.
- Power Cord 6 ft
- Power Source 120 to 240 VAC, 12 VDC (INTL fittings available); AC Adapter
- Power Status Indicator LED Indicator
- Relay Connection Spade Lugs, Male 0.25. One Relay per Channel.
- Relay Type SPST, 24 VDC, 0.5 A, N/O
- Remote Sensor Length 4
- Remote Sensor Type 1 K T/C with Mini Connector
- Response Time 30 Seconds to 63% of Full Scale
- Temperature Accuracy +/-1.8F from -22 to 122F 3.8F Remaining (+/-1C from -30 to 50C 2.1C Remaining)
- Temperature Range -100 to 0F/C; -50 to 0F/C; -50 to 50F/C; 0 to 50F/C; 0 to 100F/C; 0 to 250F/C; 0 to 500F; 50 to 100F/C
- Temperature Sensor Remote Probe
- Unit Weight 9 lbs
- Units/Pkg 1
- Warranty 12-Month Limited